




# NomoLine® Catalog




Moisture-wicking sampling lines for intubated and non-intubated patients in low- and high-humidity applications



- > Patented moisture-wicking technology allows water evaporation from the sampling line
- > Hydrophobic bacterial filter prevents contamination of the gas analyzer
- > Uniquely designed sampling line allows for usage during low-flow and small sample-size applications
- > Ergonomically designed cannulas improve patient comfort
- > Star Lumen O2 delivery tubing allows for continuous, accurate O2 supply even if kinked







## NomoLine Sampling Line Options

Division	Product Description	Part Number	Length	Image
Low Humidity or Short-term Airway Adapter Set	NomoLine LH Adult/Pediatric Airway Adapter Set, 25/box	3814	2m	
	NomoLine LH Adult/Pediatric Airway Adapter Set, 4m, 25/box	4620	4m	
High Humidity or Long-term Airway Adapter Set	NomoLine HH Adult/Pediatric Airway Adapter Set, 25/box	3827	2m	
	NomoLine HH Adult/Pediatric Airway Adapter Set, 25/box	3828	3m	
	NomoLine HH Infant/Neonatal, Airway Adapter Set, 25/box	4367	2m	
Low Humidity or Short-term Cannula with O2 Tubing	NomoLine LH Adult Nasal CO2 Cannula with O2, 25/box	3820	2m	
	NomoLine LH Adult Nasal CO2 Cannula with O2, 4m, 25/box	4624	4m	
	NomoLine LH Pediatric Nasal CO2 Cannula with O2, 25/box	3821	2m	
	NomoLine LH Adult Nasal/Oral CO2 Cannula with O2, 25/box	3824	2m	
	NomoLine LH Adult Nasal/Oral CO2 Cannula with O2, 4m, 25/box	4628	4m	
	NomoLine LH Pediatric Nasal/Oral CO2 Cannula with O2, 25/box	3825	2m	
	NomoLine LH Pediatric Nasal/Oral CO2 Cannula with O2, 4m, 25/box	4629	4m	
	NomoLine HH Adult Nasal CO2 Cannula with O2, 25/box	3833	2m	
NomoLine HH Pediatric Nasal CO2 Cannula with O2, 25/box	3834	2m		
NomoLine HH Adult Nasal/Oral CO2 Cannula with O2, 25/box	3837	2m		
NomoLine HH Pediatric Nasal/Oral CO2 Cannula with O2, 25/box	3838	2m		

Division	Product Description	Part Number	Length	Image	
Low Humidity or Short-term Cannula without O2 tubing	NomoLine LH Adult Nasal CO2 Cannula, 25/box	3817	2m		
	NomoLine LH Pediatric Nasal CO2 Cannula, 25/box	3818	2m		
	NomoLine LH Infant Nasal CO2 Cannula, 25/box	3819	2m		
	High Humidity or Long-term Cannula without O2 tubing	NomoLine LH Adult Nasal/Oral CO2 Cannula, 25/box	3822	2m	
		NomoLine LH Adult Nasal/Oral CO2 Cannula, 4m, 25/box	4626	4m	
		NomoLine LH Pediatric Nasal/Oral CO2 Cannula, 25/box	3823	2m	
		NomoLine LH Adult Single Nasal Prong CO2 Cannula, 25/box	3826	2m	
		NomoLine HH Adult Nasal CO2 Cannula, 25/box	3830	2m	
High Humidity or Long-term Cannula without O2 tubing	NomoLine HH Infant Nasal CO2 Cannula, 25/box	3832	2m		
	NomoLine HH Adult Single Nasal Prong CO2 Cannula, 25/box	3839	2m		

LH = Low Humidity    HH = High Humidity






## NomoLine-O™ Sampling Line Options

Division	Product Description	Part Number	Length	Image
Low Humidity or Short-term Airway Adapter Set	NomoLine-O LH Adult/Pediatric Airway Adapter Set, 2m, 25/Box	4454	2m	
	NomoLine-O LH Adult/Pediatric Airway Adapter Set, 4m, 25/Box	4630	4m	
High Humidity or Long-term Airway Adapter Set	NomoLine-O HH Adult/Pediatric Airway Adapter Set, 2m, 25/Box	4449	2m	
	NomoLine-O HH Infant/Neonatal Airway Adapter Set, 2m, 25/Box	4450	2m	
Low Humidity or Short-term Cannula with O2 Tubing	NomoLine-O LH Adult Nasal CO2 Cannula with O2, 2m, 25/Box	4456	2m	
	NomoLine-O LH Adult Nasal CO2 Cannula with O2, 4m, 25/Box	4634	4m	
	NomoLine-O LH Pediatric Nasal CO2 Cannula with O2, 2m, 25/Box	4457	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4459	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula with O2, 4m, 25/Box	4638	4m	
	NomoLine-O LH Pediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4460	2m	
	NomoLine-O LH Pediatric Nasal/Oral CO2 Cannula with O2, 4m, 25/Box	4639	4m	
	NomoLine-O HH Adult Nasal CO2 Cannula with O2, 2m, 25/Box	4451	2m	
NomoLine-O HH Pediatric Nasal CO2 Cannula with O2, 2m, 25/Box	4645	2m		
NomoLine-O HH Adult Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4453	2m		
NomoLine-O HH Pediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4647	2m		

Division	Product Description	Part Number	Length	Image
Low Humidity or Short-term Cannula without O2 Tubing	NomoLine-O LH Adult Nasal CO2 Cannula, 2m, 25/box	4455	2m	
	NomoLine-O LH Pediatric Nasal CO2 Cannula, 2m, 25/box	4641	2m	
	NomoLine-O LH Infant Nasal CO2 Cannula, 2m, 25/box	4662	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula, 2m, 25/Box	4458	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula, 4m, 25/Box	4636	4m	
	NomoLine-O LH Pediatric Nasal/Oral CO2 Cannula, 2m, 25/box	4642	2m	
	NomoLine-O LH Adult Single Nasal Prong CO2 Cannula, 2m, 25/box	4663	2m	
High Humidity or Long-term Cannula without O2 tubing	NomoLine-O HH Adult Nasal CO2 Cannula, 2m, 25/box	4643	2m	
	NomoLine-O HH Infant Nasal CO2 Cannula, 2m, 25/box	4664	2m	
	NomoLine-O HH Adult Single Nasal Prong CO2 Cannula, 2m, 25/box	4665	2m	

LH = Low Humidity    HH = High Humidity

## NomoLine with Luer Sampling Line Options

Division	Product Description	Part Number	Length	Image
Airway Adapter Set	NomoLine Adult/Pediatric Airway Adapter Set with Luer, 2m, 25/box	3801	2m	
Cannula with O2 Tubing	NomoLine Adult Nasal CO2 Cannula with O2 and Luer, 2m, 25/box	3807	2m	
	NomoLine Pediatric Nasal CO2 Cannula with O2 and Luer, 2m, 25/box	3808	2m	
	NomoLine Adult Nasal/Oral CO2 Cannula with O2 and Luer, 2m, 25/box	3811	2m	
	NomoLine Pediatric Nasal/Oral CO2 Cannula with O2 and Luer, 2m, 25/box	3812	2m	
Cannula without O2 Tubing	NomoLine Adult Nasal/Oral CO2 Cannula with Luer, 2m, 25/box	3809	2m	
	NomoLine Adult Single Nasal Prong CO2 Cannula with Luer, 2m, 25/box	3813	2m	

## Sampling Line Connector Versions



### HH Sampling Lines - ISA™

Long-term usage

For use with Masimo ISA and Rad-97® devices



### LH Sampling Lines - ISA

Short-term usage

For use with Masimo ISA and Rad-97 devices



### HH Sampling Lines - Medtronic® (Oridion)

Long-term usage

For use with Medtronic Microstream™ devices and modules



### LH Sampling Lines - Medtronic (Oridion)

Short-term usage

For use with Medtronic Microstream devices and modules



### Luer Sampling Lines - Water Trap

For use with luer-based devices using water traps

Not all part numbers are available for sale in Canada.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information including indications, contraindications, warnings, and precautions.

The Medtronic® trademarks are property of Medtronic and are used for identification purposes only. Use of the name, trademark, or brand does not imply any affiliation or endorsement.

**Masimo Sweden AB**  
Tel: +46 8 544 98 150  
info-sweden@masimo.com

**Masimo U.S.**  
tel 1-877-4-Masimo  
info-america@masimo.com

**Masimo International**  
tel +41-32-720-1111  
info-international@masimo.com

